





## **CBI Newsletter**

**JULY-2022** 

A Great Tool for Drug Screening Model: Monoclonal Antibody
Recognizing Wild Type Conformational Structure of p53

**p53** mutations are found in approximately half of all human cancers. **p53** is activated in response to DNA damage. Mutated p53 proteins lose their normal, tumor suppressing function.

PC14586 is a first-in-class, small molecule p53 reactivator designed to selectively bind to the crevice present in the p53 Y220C mutant protein, hence, restoring the wild-type, or normal, p53 protein structure and tumor-suppressing function. Phase 1 clinical data of PC14586, presented PMV Pharma, at ASCO demonstrated anti-tumor activity across multiple solid tumor types with a p53 Y220C mutation.


The monoclonal p53 **antibody PAb 1620** binds to a conformational epitope that disappears upon certain mutated p53 or denaturation, which provides an excellent tool for monitoring and identifying pharmacodynamic changes for this unique approach, demonstrating restoration of wildtype p53 protein structure, therefore, restores its tumor suppressor function.

#### References

- 1. Cook A, Milner J. (1990) Evidence for allosteric variants of wild-type p53, a tumor suppressor protein. Br J Cancer. 61(4):548-52.
- 2. Peter L Wang, Fiona Sait and Greg Winter (2001) The 'wildtype' conformation of p53: epitope mapping using hybrid proteins. Oncogene 20 (18): 2318-24
- 3. https://ir.pmvpharma.com/news-releases/news-release

Figure. MCF-7 breast cancer cells were treated with (D-I) or without (A-C) 50 nM dox for 24 hrs, incubated with IgG2a isotype control or anti-p53 PAb-1620 diluted 1:1000. Donkeyanti-mouse iFluor<sup>TM</sup> 488 (green, B, E, H) was used in a dilution of 1:1000, and DAPI was applied for detection of nuclei (blue, A, D, G).

### Representative data





# Our available products:

PRODUCT	CLONE	CATALOG#
Anti-p53 Unconjugated	PAB122	102101
Anti-p53 WT Unconjugated	PAB1620	102201

### **JULY-2022**

### **NEW PRODUCTS**

Product	Size	Clone	Catalog #
CD10 APC	100Test	C-1A12	114945
CD45RO PerCP-Cyanine5.5	100Test	UCHL1	107865
IgG1 PE-CF594	100Test	MOPC-31C	1023205
CD34 PE-CF594	100Test	4H11	1089205
CD23 APC	100Test	HD-50	109545
CD27 FITC	100Test	C-813.21	125715
CD27 APC-iFluor700	100Test	C-813.21	122025
CD27 PE	100Test	C-813.21	125725
MPO APC	100Test	8/9	112345

### 200 tests & ASR & free test size available!













Caprico Biotechnologies is committed to



9 400 Pinnacle Way, Suite 430 Norcross, GA 30071

100% customer satisfaction.

